

QUICK & EARLY detection of COVID-19 Infection



RESZON COVID-19 RAPID ANTIGEN TEST (SALIVA SAMPLE)

CE IVD FOR PROFESSIONAL USE ONLY, NOT FOR HOME TESTING



RESZON COVID-19

RAPID ANTIGEN TEST (SALIVA SAMPLE)

QUICK & EARLY detection of COVID-19 infection


Test Description

This is an immunochromatographic assay designed for the qualitative detection of specific SARS-CoV-2 antigen in human saliva specimen. It is intended to be used by the professionals as a screening test and as an aid in the early diagnosis of SARS-CoV-2 infection.


Features

Platform	Immunochromatographic Assay
Format	Cassette
Detection	Specific SARS-CoV-2 antigen
Specimen	Saliva
Assay Time	15-20 minutes
Shelf Life	18 months

Procedures (Simplified Version)



Result Interpretation



Ordering Information

Cat. No	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp	Packing Size
RCV-RD0201S-10	COVID-19 Rapid Antigen Test (with Saliva Collection Set)	Saliva sample	Cassette	Coronavirus disease COVID-19	4-30°C	10 T/ Kit

COVID-19

